Product Specification and Description Sheet

Product name 5-Bromo-5-nitro-1,3-dioxane

Code #BND-D

Purity (GC) ≥99,5%

Formula $C_4H_6BrNO_4$

CAS# 30007-47-7

Molecular weight 212.0

Appearance white to slightly yellowish dry free

flowing crystalline powder

Melting point 58,3-61,1°C

Boiling point 210°C

Absorbance maximum (SP) 224,5 nm (solution in EtOH)

Solubility in water 0.5% in H₂O >18M Ω (at 29°C in 2 hours)

Rest H_2O (Karl Fischer titration) $\leq 0.5\%$

Storage stability >>2 years at 15 - 25°C (in tightly closed

glass or HD-PE/PP/PET)

Activity very strong antimicrobial preservative

highly effective against bacteria, yeast and


fungi

Specific conc. effective in test with

Gram-positive bacteria: Bacillus subtilis ATCC 6633, Streptococcus

faecalis (clinical strain)

100 ppm (0,01%)



SDT GmbH • Pascalstraße 17 52499 Baesweiler • Germany www.sdt-reagents.de

Gram-negative bacteria: Esherichia coli ATCC 11229 (extremely resistant strain), Pseudomonas aeruginosa ATCC 14207, Enterobacter clocae ATCC 23355

300 ppm (0,03%)

Yeast/Fungi: Candida albicans ATCC 10231, Aspergillus niger ATCC 6275, Penicillium funiculosum ATCC 9644, Chaetomium globosum ATCC 6205 extremely resistant strain)

100 ppm (0,01%)

Suggested mechanism of action

presumably oxidizes thiol groups in vitally important microbial enzyme systems

Recommended concentrations

0,02-0,03% (simple salt buffers) – to - 0,05% (buffers containing low to moderate concentrations of proteins, sugars or detergents) – to - 0,1-0,12% (products containing high protein concentrations, liquid protein concentrates, ready-to-use IVD reagents)

pН

5,0-6,0 in pure water has no detectable influence on pH does not alter pH in buffers through 3,5-9,5


Applicable pH range

3,5-9,5

Compatibility and Interferences

compatible with most common liquid diagnostic formulations; does not interfere with antigen-antibody reactions, most enzyme assays, immunoassays, immunosensors, restriction enzymes, PCR,

electrophoresis, nucleic acid and protein purification protocols.



SDT GmbH • Pascalstraße 17 52499 Baesweiler • Germany www.sdt-reagents.de

interferes with cysteine and other free thiol

reagents

Restrictions not stable at temperatures >45°C

slightly corrosive to metals

strong reducing agents may lower biocide

potential

Safety/Toxicity harmless at recommended concentrations;

not carcinogenic or mutagenic (Ames-test

neg.)