

Specification Sheet

Product:	Drying Buffer/Stabilizer, biotin-free BSA based, ready-to-use
Code:	#DBS-B
Protein concentration :	100 mg/ml
Buffer base:	0,15M PBS pH 7,3
Preservative:	BND, 550 ppm
Lot:	538892
Storage:	+10°C to +25°C in darkness (recommended) +2°C to +8°C (alternative storage regimen giving longer 2-year stability; prior to use warm up at RT/+37°C to allow for dissolving precipitate, if the latter has been formed in storage)
Expiration date:	12/2007+

Recommended use: ready-to-use reagent for preparing dried non porous and porous immunosorbents that shall be stable real-time at room/ambient temperature